ATOS BLOCK

MMTech LTDA

Rua Doutor Procópio de Toledo Malta - 62
Morada dos Deuses - CEP 13562-291

São Carlos-SP – Brasil

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Product Application

The Dental Surgeon will mold or scan the patient's oral cavity, design the part in dental CAD software and mill the block according to the desired restorative treatment.

The restorations produced must be placed with an adhesive composite luting agent in combination with a dentine-enamel bonding agent. Available in size 14 L. Universal holder (fitting) - Cerec Sirona® type.

Core and cavity preparation:

As a general rule, the core/cavity preparation should be performed in accordance with the conventional rules for all-ceramic restorations. That means rounding off interior corners and edges and using a shoulder preparation with rounded interior edges or a chamfer preparation. The temporary restoration and any contamination should be mechanically removed (e.g.: using pumice and a polishing brush) and the area subsequently rinsed with water. Observe the minimum thicknesses stated below. Note: Do not create any

eugenol based materials as they may inhibit proper cement/adhesive curing.

Product description

ATOS BLOCK is a CAD/CAM nano-hybrid composite block for permanent single-tooth restorations. It is recommended for restorations of all classes in anterior and posterior teeth. The composite has fluorescence, radiopacity and contains 77% to 86 % by weight of inorganic filler particles in a highly polymerized, Bisphenol A Free [BPA Free] polymer matrix. ATOS BLOCK has a high degree of conversion, significantly higher than conventional light-curing composite resins. These factors allied to adequate particle size distribution result in a composite with high hardness, mechanical strength, and wear resistance, allowing parts to be obtained with excellent finishing, smooth and high-gloss surfaces.

Indications:

- Unitary full crown
- Unitary full crown cemented over abutment
- Inlays and Onlays
- Veneers

Shade selection:

Use the VITA® shade system, for guide, to determine the shade against the cleaned but unprepared tooth preferably in daylight conditions.

Shades (two degrees of translucency):

High Translucent (HT): A1, A2, A3, B1 Low Translucent (LT): A1, A2, A3, A3.5, B1, C2, B3

Product Application

The Dental Surgeon will mold or scan the patient's oral cavity, design the part in dental CAD software and mill the block according to the desired restorative treatment.

The restorations produced must be placed with an adhesive composite luting agent in combination with a dentine-enamel bonding agent. Available in size 14 L. Universal holder (fitting) - Cerec Sirona® type.

Core and cavity preparation:

As a general rule, the core/cavity preparation should be performed in accordance with the conventional rules for all-ceramic restorations. That means rounding off interior corners and edges and using a shoulder preparation with rounded interior edges or a chamfer preparation. The temporary restoration and any contamination should be mechanically removed (e.g.: using pumice and a polishing brush) and the area subsequently rinsed with water. Observe the minimum thicknesses stated below. Note: Do not create any

undercuts.

Crowns, inlays and onlays (see Fig. 1- 3):

- Preparation angle of 4°- 6°
- Minimum occlusal wall thickness 1,5 mm
- Minimum buccal wall thickness 0,8 mm
- Minimum wall thickness below load-bearing cusp 1,5 mm
- Minimum cervical wall thickness 0,8 mm

Veneers (see Fig. 4):

- Minimum cervical wall thickness 0,4 mm
- Minimum labial wall thickness 0,5 mm
- Minimum incisal wall thickness 0,5 mm

Cement space:

70 µm (± 10 µm)

Milling process:

Select the correct block size for the planned restoration as well as the milling parameters for ATOS BLOCK or for composites. When doing so, pay attention to the software settings of the respective CAD/CAM systems and user manuals from the CAD/CAM milling machine manufacturer.

Check the processed restoration for defects such as fissures or material chipping. If the restoration displays any defects, it must be rejected.

Note: ATOS BLOCK is an already completely cured composite and must not be fired.

Preparing the restoration:

For an optimal bond, blast the luting surface of the restoration with aluminium oxide [25 µm - 50 µm] at 1,5 bar - 2 bar. Use a suction de-vice to remove the dust produced.

NOTE: Do not etch with hydrofluoric acid or phosphoric acid.

Clean the pretreated restoration using an ultrasonic bath or steam cleaner and then dry with oil-free air. Final cleaning with medical alcohol is possible. Treat the luting surface with a silane. Apply adhesive (that contain MDP) to the inside of the work and allow it to work for 60 seconds. Then dry for 5 seconds with oil-free air. The instructions for use must be observed and followed.

Adhesive luting of the restoration:

The processed restoration must be luted with a composite-based, adhesive luting agent. The instructions for use must be observed and followed. Alternatively, it is also possible to use a purely light- cured, flowable composite if the restoration has thin walls [≤ 2 mm].

Preparation of the trunnion on implants:

The post must receive dental adhesive that

contains MDP throughout its cementable area, wait 60 s and light cure for 15 s.

Occlusal adjustments on the workpiece:

Once the part is ready, it should be tried in mouth before cementation. Any necessary occlusal adjustments can be made with fine-grain diamond drills.

Polishing:

The restoration can be finished and polished with standard compo- site polishers (e.g. medium, fine or extra fine diamond grinding tools, polishers), this can be done either intraorally or extraorally.

Modifications:

ATOS BLOCK is a nanohybrid composite CAD/CAM restorative material which can be modified, characterized, or repaired at any time using a composite (e.g., Atos). Roughen the surface of the resto- ration (for example with a diamond bur) and apply a suitable adhesive system (that contain MDP) in accordance with the instructions for use. Then apply the composite (e.g., Atos) in accordance with the respective instructions for use.

Presentation forms

Printed on the product packaging. Pack containing 5 units of blocks. Pack containing 1 to 20 units of blocks.

NOTE: Instructions for use for the professional can be found on the Smart Dent website. Use the QR Code on the packaging.

Basic composition:

Contains methacrylate monomers, stabilizers, initiator, fillers, silane and pigments.

Expiration date:

Printed on the product packaging.

Storage conditions:

The product must remain in its original packaging. Keep the product in a cool place. Protect from direct sunlight. Store the product at temperatures from 5 °C to 27 °C / 41 °F to 80 °F. Do not freeze the product.

Handling Conditions

- The product is indicated for professional use only.
- Read the instructions and instructions for use carefully.
- Do not use the product if it is out of date.

- To dispose of the product, follow the legislation of your country.
- Clean instruments with alcohol after use.
- The use of gloves is recommended.

Warnings

Do not use the product if it is out of date. Keep out of the reach of children. Check the compatibility of Atos Block milling with your equipment manufacturer.

Precautions

Professional use only

- This product is indicated only according to its instructions for use.
- Avoid the use of eugenol based lining or provisional materials, as these interfere with the polymerization of the material.
- Avoid use in patients with a history of allergy to composite resins based on methacrylates. In case of allergy, discontinue use and refer the patient for medical assistance.
- When using a light curing device, eye protection must be used.
- Light curing equipment must have its power checked with a radiometer frequently.

Contraindications and adverse effects

The product is contraindicated for patients allergic to methacrylates. The product contains methacrylate monomers that can cause mild allergic or irritant reactions in patients sensitive to these substances. In case of allergic reactions to the product, use must be suspended.

This material has been manufactured for dental used only and must be handled according to the use instructions and manufacturer's manual for CAD/CAM milling machines. The manufacturer is not responsible for damage caused by other uses or incorrect handling and equipment(s). In addition, the user is obliged to check, before use and under his responsibility, if this material is compatible with the intended use, especially when this use is not indicated in the use instructions or in the milling machine manual. Data descriptions do not constitute any type of guarantee and, therefore, have no link.

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