# PERFIL TÉCNICO

ISO 4049

# ATOS

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Indústria Brasileira - Made in Brazil

Marca SMART DENT® - Rev: 00

 SMART DENT

SAC +55 16 33743950 | [www.smartdent.com.br](http://www.smartdent.com.br)  
[sacdentistica@smartdent.com.br](mailto:sacdentistica@smartdent.com.br)

# TECHNICAL PROFILE

ISO 4049

# ATOS

**Nano-hybrid composite resin with Nano Zirconium Silicate**

 SMART DENT

# SUMMARY

<b>Introduction</b>	2
<b>Main Technological Features</b>	3
Nano-hybrid	3
Polishing Shine	3
Photopolymerization system	3
Colors and opacity	4
Fluorescence, opalescence and radiopacity	4
Monomers and structuring agents	5
Structural resistance	5
Composition	5
Indication of Use	6
<b>Physicochemical and Mechanical Properties</b>	7
Flexural strength	7
Surface roughness	8
Microhardness	8
Compressive strength	9
Fracture toughness	10
Conversion degree	10
Modulus of elasticity	11
<b>Specialist Report Clinical</b>	12
<b>Cases</b>	13
Aesthetic and functional reconstruction of a dentally treated tooth - Prof. Dr. Maciel Jr.	13
Mimicking teeth naturally - Prof. Dr. Maciel Jr.	14
Aesthetic properties of Atos in a clinical case of diastema closure - Prof. Dr. Victor Feitosa High-	15
standard aesthetic and functional resin restoration - Prof. Dr. Renato Voss Rosa	16

# INTRODUCTION

Atos is a light-curing nano-hybrid composite resin, indicated for restorations of all classes in anterior and posterior teeth using the incremental system.

This universal dental restoration composite has fluorescence, opalescence and radiopacity. The formulation of the Atos monomeric matrix was based on the main structuring methacrylic monomers. Its inorganic composition is composed of a distribution of charges in which Zirconium Nano Silicate stands out. The total content of inorganic charge by weight is  $79.0 \pm 2\%$  and  $60.0 \pm 2\%$  by volume. The average particle size is 0.3 microns.

Atos has a photo system activated by blue light at a wavelength of 400 to 500 nm for Halogen or LED light and with power from  $600 \text{ mW/cm}^2$  to promote its polymerization.

Atos restorative composite is available in 20 shades. Its range of colors and opacity variations were developed and based on the highest mimicry technology for dental aesthetics, allowing success in the layering technique, providing maximum aesthetic performance and the ability to mimic natural teeth.

Atos resin therefore provides high aesthetics for natural teeth, adjusts the harmony of the smile, as well as restores masticatory functions and contributes to the well-being of dental patients.


# MAIN TECHNOLOGICAL FEATURES

TECHNICAL PROFILE - ATOS BLOCK

SMART DENT

## NANO-HYBRID

Its inorganic composition is composed of silanized fillers manufactured in Germany and has high performance, such as Barium-Aluminum glass in the range of  $0.7\ \mu\text{m}$ , Silica nanoparticles in the range of  $35\ \text{nm}$  and mainly Zirconium Nano Silicate in the range of  $100\ \text{nm}$ , which maximize wear resistance. The uniform distribution of particles from  $35\ \text{nm}$  to  $1.0\ \mu\text{m}$  contributes to high mechanical properties, excellent handling, high aesthetics and especially shine when polished.

## POLISHING SHINE

The spheroidal shape, the size of 100 nanometers and the added percentage of the Zirconium Nano Silicate filler, bring to the restorative composite the power of high brightness and surface smoothness that is visually close to that of a feldspathic ceramic and with prolonged polishing retention.

Scanning Electron Microscopy (SEM) photomicrograph of the spheroidal particles of the ATOS restorative composite. Magnification of 10,000 times.

Image kindly provided by PUCPR

## PHOTOPOLYMERIZATION SYSTEM

The ATOS blue light polymerization system provides professionals with greater safety regarding the exposure time of the photopolymerizer, especially in relation to its main competitors, for example: enamel colors can be polymerized in 20 to 30 seconds and dentin in 30 to 40 seconds, as they already meet the ISO requirements.

4049. 'NOTE: the longer time specified gives the formed polymer greater properties: In this way, it is possible to polymerize the material after working on its adequate positioning and shape of the restoration

## COLORS AND OPACITIES:

The VITA® Classical shade scale can be used as a reference for choosing traditional shades. Dentin shades exhibit greater opacity and a more saturated chroma, while enamel shades have greater translucency, adjusted to diffuse the dentin substrate. The effect shades and shades for whitened teeth were developed based on recommendations from dental surgeons specializing in aesthetics. All shades are radiopaque except for Translucent OPL, due to its translucency and high optical effect.

Atos has 20 color options that meet the needs of the simplest to the most complex aesthetic restorations.

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
## FLUORESCENCE, RADIOPACITY AND OPALESCENCE:

**Fluorescence** is reproduced by special pigments, which absorb energy from UV light and emit it through visible light. The Atos system is balanced to reproduce this effect, masterfully imitating the fluorescence of natural teeth.

The **radiopacity** of the Atos system, in addition to being in accordance with ISO 4049, provides the dentist with a safe and accurate diagnosis.

**Opalescence** is a unique characteristic of natural tooth enamel that absorbs blue colors and allows orange colors to pass through. This optical effect is found in Atos enamel colors, although in a subtle way, but it is extremely evident in the Translucent OPL (Opalescent) color.

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
## STRUCTURING MONOMERS

The organic matrix of this resin is formed by the main structuring monomers in the development of composites in the world, manufactured in the USA and with high molecular weight. They are Bis-Phenol A di-Glycidyl Methacrylate (BisGMA), 1, 12-dodecanediol dimethacrylate (DDDMA), Bis-Phenol A di-Glycidyl Methacrylate ethoxylated (BisEMA), UDMA (Urethane dimethacrylate) and with a slight diluent part of Triethylene glycol dimethacrylate (TEGDMA). These monomers perform important functions, such as reduced polymerization shrinkage factor, color stability, and increased service life. thixotropy and ease of handling.

## MATERIAL STRENGTH

The combination of nanometric particles with micrometric particles in the formulation reduces the space between the filler particles in the monomeric solution, increasing the percentage of filler which, inserted into the structuring monomers, considerably increases the physical and mechanical properties of the resin, providing high resistance, long-lasting restorations with the best possible aesthetic results.

Scanning Electron Microscopy (SEM) photomicrograph of the resin particle combination. Atos presents uniformity and does not present discrepant (large) particles. **Magnification of 10,000 times.** Image kindly provided by FOP-Unicamp

## COMPOSITION

Atos resin is composed of a matrix with methacrylate monomers of BisGMA, BisEMA, UDMA, DDDMA and a small portion of TEGDMA. Photoinitiator system, silane and stabilizers. The inorganic filler is composed of Barium-Aluminum glass, Silica nanoparticles and mainly Zirconium Silicate and pigments.

### Cores

#### Enamel:

EA1, EA2, EA3, EA4, EB1, EB2, WE, YOURSELF

#### Dentine:

DA1, DA2, DA3, DA3.5, DA4, DB1, DB2, WD, XWD

#### Effect:

Opaque, Transparent Clear, Opalescent OPL

### Light curing time

Cor	Thickness	Time
Dentine	2 mm	30 - 40 s
Enamel	2 mm	20 - 30 s
Translucent	2 mm	20 s
Opaque	2 mm	20 s

## INDICATIONS FOR USE

- Restorations of permanent and deciduous teeth, anterior and posterior, of classes I, II, III, IV and V (all classes)
- Direct composite resin veneers
- Bonding of tooth fragments
- Anchoring of teeth (splinting)
- Correction of tooth shape (e.g. conoids)
- Reduction or closure of diastemas
- Production of filling cores
- Non-carious cervical lesions
- Porcelain and/or composite repairs
- Indirect restorations, such as inlays, onlays and veneers
- Cementation of prosthetic parts (according to the thermo-modified resin technique)

## WORD FROM THOSE WHO KNOW:


The Atos restorative composite was tested at our Institute and is equivalent to the main international brands, it has high-level anatomical restorations with excellent handling...

**An excellent national option.**

**PROF. DR. JOSÉ MONDELLI**

PROFESSOR OF MONDELLI INSTITUTE OF ONDOTOLOGY - BAURU-SP

Laboratory tests carried out to date show excellent results from Atos resin, a national technology showing good performance!!

**PROF. DR. BRUNO REIS** – UBERLÂNDIA-MG

# PHYSICOCHEMICAL AND MECHANICAL PROPERTIES


Tested by **FACPP** and **FOP-UNICAMP** universities

## FLEXIBLE STRENGTH

Flexural Strength is one of the most important mechanical tests for dental composites as it evaluates the fracture resistance of the restorative material as well as its intrinsic properties (monomeric matrix, inorganic fillers and photoinitiators) after photopolymerization, in relation to different masticatory forces.

The specimens of each resin were prepared in a bipartite stainless steel matrix with specific dimensions (25x2x2 mm) in the shape of a rod. After three photoactivations of 20 s (one on each side and one in the center), with irradiance of 1200 mW/cm<sup>2</sup>, the resins were subjected to a three-point bending test in an EMIC 23-2S universal testing machine (EMIC-Instron, São José dos Pinhais, Brazil) until fracture, according to ISO 4049.

**Figure 1:** Flexural strength at 3 points (mean and SD in MPa) in 4 different composites (n=10). Statistical analysis ANOVA and Tukey test; p<0.05.


**Source:** Rifane TO, Cordeiro KE, Silvestre FA, Alves AHC, Nobre CFA, Araújo-Neto VG, Giannini M, Feitosa VP. Paulo Picanço College of Fortaleza (FACPP) and Faculty of Dentistry of Piracicaba (FOP-Unicamp) 2021.

**ATOS** presented surprising **FLEXIBLE STRENGTH** values, being **SUPERIOR TO ALL** commercial **BRANDS** of competing resins tested.

## SURFACE ROUGHNESS

The longevity of aesthetic and functional restorations depends directly on their resistance to abrasion. Analysis of the roughness test before and after brushing reveals whether the restorative composite has appropriate wear resistance properties. High-quality resins should have high wear resistance and therefore low roughness after brushing.

The surface roughness of each specimen was evaluated using a roughness meter (Hmv-2 Shimadzu, Japan), with a cut-off of 0.08 mm and speed of 0.1 mm/s, before (initial) and after brushing for a total of 100,000 cycles on a brushing and abrasiveness machine 4 plus (Odeme, Luzerna, Brazil). The reading considered was the arithmetic mean between peaks and valleys (Ra), with three readings being taken on each specimen in different directions, according to ISO 6872.

**Figure 2:** Surface roughness (mean and SD in Ra) in 4 different composites (n=5), after polishing with 2000 Mesh sandpaper. Statistical analysis ANOVA and Tukey test;  $p < 0.05$ .


**Source:** Rifane TO, Cordeiro KE, Silvestre FA, Alves AHC, Nobre CFA, Araújo-Neto VG, Giannini M, Feitosa VP. Paulo Picanço College of Fortaleza (FACPP) and Faculty of Dentistry of Piracicaba (FOP-Unicamp) 2021.

The **ATOS**, Z250XT and Charisma resins obtained the **BEST RESULTS**, presenting lower surface roughness after polishing with sandpaper.

## MICROHARDNESS

Hardness is a property that expresses the mechanical strength and wear resistance of the composite resin. This property is influenced by the degree of monomeric conversion after photoactivation. It is also an indirect way of evaluating gloss maintenance. The higher the hardness, the better its physical, chemical and mechanical properties.

For microhardness analysis, a SE1700 microhardness tester (Kosak lab, Akita, Japan) with a Knoop type indenter and a static load of 50 grams for 10 seconds was used. Three indentations were made in triplicate, separated by a distance of 100 m. These were made on the upper surface of each test specimen.

**Figure 3:** Microhardness (mean and SD of Knoop Hardness number) in 4 different composites (n=5) in triplicate. Statistical analysis ANOVA and Tukey test; p<0.05.


**Source:** Rifane TO, Cordeiro KE, Silvestre FA, Alves AHC, Nobre CFA, Araújo-Neto VG, Giannini M, Feitosa VP. Paulo Picanço College of Fortaleza (FACPP) and Faculty of Dentistry of Piracicaba (FOP-Unicemp) 2021.

The composite resins **ATOS** and **Z250XT** presented statistically **higher hardness values** than the **other brands** tested.

## COMPRESSION STRENGTH

The mechanical property of compressive strength demonstrates the ability of the composite resin to withstand vertical stresses, as the forces transmitted to the restorations during the chewing act can trigger fractures of the material.

The higher the compressive strength value, the longer the composite will last. Cylindrical samples were prepared for each material using a silicone rubber mold with dimensions of 4 mm (in diameter) × 8 mm (in height). The mechanical compression test was performed on the universal testing machine.

**Figure 4:** Compressive strength (mean and SD in MPa) in different composites (n=5). Statistical analysis ANOVA and Tukey test; p<0.05.


**Source:** Rifane TO, Cordeiro KE, Silvestre FA, Alves AHC, Nobre CFA, Araújo-Neto VG, Giannini M, Feitosa VP. Paulo Picanço College of Fortaleza (FACPP) and Faculty of Dentistry of Piracicaba (FOP-Unicemp) 2021.

**ATOS** presents **EXCELLENT COMPRESSION RESISTANCE**, comparable to well-known products on the market and superior to Forma resin.

## FRACTURE TENACITY

Fracture toughness or crack propagation resistance is a measure that describes the ability of a material, with a crack or structural defect, to resist fracture. This measure, therefore, indicates the amount of force required to cause the extension of the crack in the material. The higher the value of the resistance, the less likely the crack will propagate in a material.

The specimens of each resin were made in a bipartite stainless steel matrix with specific dimensions (25x2x2 mm) with a central slit of 2.58 mm. After 3 photoactivations of 20 s on each side, with irradiance of 1200 mW/cm<sup>2</sup>, according to ASTM E-399 Standard, the resins were tested in the universal testing machine and subjected to compression load in the center of each specimen in the 3-point test until fracture.

**Figure 5:** Fracture toughness at 3 points (mean and SD in MPa/m<sup>2</sup> in 4 different composites (n=10). Statistical analysis ANOVA and Tukey test; p<0.05.


**Source:** Rifane TO, Cordeiro KE, Silvestre FA, Alves AHC, Nobre CFA, Araújo-Neto VG, Giannini M, Feitosa VP. Paulo Picanço College of Fortaleza (FACPP) and Faculty of Dentistry of Piracicaba (FOP-Unicamp) 2021.

**ATOS** showed balance in its composition, and it was observed that only in this study all the results are similar. ATOS is always **EQUAL OR BETTER** in the other tests.

## DEGREE OF CONVERSION

The degree of conversion of a composite resin is given by how much monomer was converted into polymer, after the polymerization process of the dental composite.

The degree of conversion was analyzed by a Spectrum Frontier FTIR spectrophotometer (Perkin-Elmer Corp., Norwalk, United States) after 20 s polymerization of the composite resins, equipped with an attenuated total reflectance (ATR) device. The peaks at 1609 cm<sup>-1</sup> (aromatic bonds) and 1639 cm<sup>-1</sup> (aliphatic double bonds C=C) were verified. This aromatic band was used as an internal standard. The ratio between peaks was calculated for both the polymerized and unpolymerized material.

**Figure 6:** Degree of conversion (mean and SD in %) in 4 different composites (n=5). Statistical analysis ANOVA and Tukey test; p<0.05.


**Source:** Rifane TO, Cordeiro KE, Silvestre FA, Alves AHC, Nobre CFA, Araújo-Neto VG, Giannini M, Feitosa VP. Paulo Picanço College of Fortaleza (FACPP) and Faculty of Dentistry of Piracicaba (FOP-Unicamp) 2021.

the **ATOS** composite presented a **SIMILAR DEGREE OF CONVERSION** to the other products tested and is within an excellent conversion range. The **EXCELLENT RESULTS** in the other tests are evidenced.

## ELASTICITY MODULE

The term modulus of elasticity or flexural modulus describes the rigidity of the resin and is important for determining the resistance to occlusal forces and wear resistance, which reflect an intrinsic property of the material to resist pressures imposed by masticatory forces. The higher the elastic modulus, the greater the resistance of the composite to fracture.

The specimens of each resin were made in a bipartite stainless steel matrix with specific dimensions (25x2x2 mm] in the shape of a rod. After 3 photoactivations of 20 s (one on each side and one in the center), with irradiance of 1200 mW/cm<sup>2</sup>, the resins were subjected to a three-point bending test in an EMIC 23-2S universal testing machine (EMIC-Instron, São José dos Pinhais, Brazil) until fracture, according to ISO 4049.

**Figure 7:** 3-point flexural modulus (mean and SD in GPa) in 4 different composites (n=10). Statistical analysis ANOVA and Tukey test; p<0.05.

**Source:** Rifane TO, Cordeiro KE, Silvestre FA, Alves AHC, Nobre CFA, Araújo-Neto VG, Giannini M, Feitosa VP. Paulo Picanço College of Fortaleza (FACPP) and Faculty of Dentistry of Piracicaba (FOP-Unicamp) 2021.

**ATOS** presented an elastic modulus **COMPARIABLE TO OR SUPERIOR** to international products already known on the market.


# EXPERT REPORTS

Excellent handling capacity and good opacity between colors. Atos resin appears as an excellent option for the day-to-day use of the **most experienced clinician and the beginner as well.**

PROF. DR. RAPHAEL MONTE ALTO - NITERÓI-RJ

Atos resin has excellent sculptability, making it easy to obtain adequate anatomy. As a professor of functional dental anatomy, **I recommend Atos for anterior and posterior restorations.**


PROFª. DRª. TININHA GOMES - JOINVILLE-SC

Seeing a **national material with this quality is really cool!** I loved a lot of things, but what caught my attention was the opacity of the material... just right!

PROF. DR. THIAGO OTTONI - BLUMENAU-SC

What impressed me most about ATOS was its **excellent consistency, allowing for agility in the anatomization of restorations.**

PROFª. DRª. SIMONE ALBERTON - PROFESSORA DO CEOM - PASSO FUNDO-RS

# CLINICAL CASE

## AESTHETIC AND FUNCTIONAL RECONSTRUCTION OF ENDODONTIC TREATED TOOTH

Author: Prof. Dr. Maciel Jr.

1. Initial case

2. Preparation

3. Reconstruction of nucleus and near face

4. Photopolymerization

5. Immediate End

6. Clinical case completed - polished and hydrated

# MIMICKING TEETH NATURALLY

Author: Prof. Dr. Maciel Jr.

1. Initial case

2. Preparation and selective acid etching

3. Application of the adhesive system

4. Immediate end

**PROF. DR. MACIEL JR**

PROFESSOR FAEPO UNESP - ARARAQUARA-SP

The **Atos** light-curing composite resin is **excellent**, well above the world standard average!

# AESTHETIC PROPERTIES OF ATOS IN A CLINICAL CASE OF DIASTEMA CLOSURE

Author: Prof. Dr. Victor Feitosa

1. Choosing the color of the Atos resin

2. Initial case, patient with diastema - Lateral view

3. Initial case, patient with diastema - Front view

4. Application of phosphoric acid - Teeth 13 to 23

5. Clinical case performed - Side view

6. Clinical case performed - Front view

“Atos composite resin has **excellent physical and chemical properties, meeting international standards.** In addition, it stands out for its excellent handling, facilitating the achievement of restorations with excellent aesthetics, shine and polishing.”

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# RESTORATION WITH HIGH AESTHETIC AND FUNCTIONAL STANDARD RESIN


Author: Prof.Dr. Renato Voss Rosa

1. Initial case, amalgam replacement

2. Tooth preparation

3. Application of acid attack


4. Beginning of restoration after application of adhesive

5. Restored tooth

6. Tooth finished with Atos WD simulating natural tooth roughness

PROF. DR. RENATO VOSS ROSA - CURITIBA-PR

A restorative material developed and produced in Brazil is encouraging news. **Laboratory tests suggest high performance**, which is promising and has generated good expectations.

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