PERFIL TÉCNICO

ISO 4049

3em1 BULK FILL
INCREMENTAL
UNICROMÁTICA
Bulk fill | Incremental | Unichromatic

ATOS *Unichroma*

Resina composta fotopolimerizável unicromática bulk fill

 SMART DENT

SUMÁRIO

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Indústria Brasileira - Made in Brazil

Marca SMART DENT® - Rev: 00

 SMART DENT

SAC +55 16 33743950 | www.smartdent.com.br
sacdentistica@smartdent.com.br

TECHNICAL PROFILE

ISO 4049

3ⁱⁿ1 UNICROMATIC
INCREMENTAL
BULK FILL

Bulk fill | Incremental | Unichromatic

ATOS *Unichroma*

Bulk fill unichromatic light-curing composite resin

 SMART DENT

SUMMARY

Introduction	3
Description	4
Main Technological Features	5
Mimicry	5
Dentistry and Pediatric Dentistry	6
Photopolymerization System	6
Nano-hybrid	6
Fluorescence, opalescence, and radiopacity	6
Structuring Monomers	7
Material Strengths	7
Composition	7
Monochromatic	7
Indications for Use	8
Advantages	8
Physicochemical and Mechanical Properties	9
Depth of polymerization	9
Degree of conversion	10
Flexural strength	10
Compressive strength	11
Specialist Report	12
Clinical Cases	13
Restoration with Atos Unichroma composite resin of fractured central incisor - Prof. Dr. Fabio Sene Atos Unichroma	13
bulk fill - simplicity, efficiency and aesthetics in posterior teeth - Prof. Dr. Carlos Francci	14
Also discover the Smart Dent aesthetics line	15

INTRODUCTION

Color is an important element for clinical success, especially for aesthetic restorations. Color difference in restorations made on natural teeth with composite resins is one of the most common reasons that affects clinical appearance and patient dissatisfaction, making restoration replacement a frequent phenomenon in clinical practice.

Due to the variety of colors found in individuals, companies in the dental sector have long created composites with different types of colors, opacities, optical effects, etc. to try to meet the need to restore different shades of teeth. However, the color matching, the choice of tone and opacity, and the restoration with different colors considerably increase the time in the chair, which is expensive for dentists and uncomfortable for patients.

To make dentists' lives easier, SMART DENT has created the ATOS UNICHROMA light-curing composite resin, which, with a single shade and opacity, covers all the shades of the VITA® Classical scale (from A1 to D4) as well as whitening teeth. Choosing shades, learning about each shade of each resin from each manufacturer, buying different brands of resins... is no longer necessary!

ATOS UNICHROMA, with just one color, has the ability to combine all colors and opacities thanks to SMART CARE PHOTOCROMA technology. The technology is obtained by controlling the refractive index between the structuring monomers and high-performance fillers, the intelligent photoinitiator system and the mirroring of the light contained in the composite.

ATOS UNICHROMA: a resin for all cases and colors!

DESCRIPTION

THE WORLD'S FIRST COMPOSITE THAT IS BULK, INCREMENTAL AND MAINLY UNICROMATIC (SINGLE COLOR FOR ALL CASES FOR ALL COLORS)! IT'S BRAZILIAN! IT'S SMART!

Atos Unichroma composite resin is a light-curing dental restorative composite with a single color and opacity for all dental shades, called unichromatic resin (universal color); indicated for restorations of all classes in anterior and posterior teeth using the incremental or bulk fill system. The composite is nano-hybrid and has fluorescence, opalescence and radiopacity. Its inorganic filler is composed of Barium-Aluminum glass, Silica nanoparticles and Zirconium Nano Silicate. The average particle size is 0.3 microns, with a variation in the size of silanized particles from 35 nm to 0.7 μm , characteristics that contribute to high mechanical properties and high aesthetics, especially polishing shine. The total inorganic filler content by weight for Atos Unichroma resin is 74% to 80% (56% to 61% by volume).

Atos Unichroma resin is composed of a matrix with structuring dimethacrylate monomers of BisGMA, BisEMA, UDMA, DDDMA and TEGDMA. Its polymerization system is through photoactivation with blue halogen light or LED. In this way, the material is polymerized after working on the appropriate positioning of the composite resin and the shape of the restoration. It has low polymerization shrinkage. No covering resin is required! This material was developed and based on the highest mimicry technology, allowing success with just a single resin color in a true chameleon effect. If there are surrounding walls, Atos Unichroma mimics the color of the tooth, due to its excellent refractive index with color scattering, it is possible to copy the dental substrate. The composite with the aforementioned characteristics provides aesthetic performance, mechanical resistance and color mimicry capacity. Atos Unichroma resin therefore provides aesthetics to natural teeth as well as restoring chewing functions, ease of use and agility in restorations and contributes to the well-being of dental patients.

MAIN TECHNOLOGICAL FEATURES


TECHNICAL PROFILE - ATOS UNICHROMA

SMART DENT

MIMICRY:

Atos Unichroma is a universal color resin. Yes! A single color and opacity for all cases and classes. Whether incremental or Bulk Fill with up to 5 mm. Excellent for the Dentistry and Pediatric Dentistry areas.

The restorative composite has in its composition, an exclusive SMART CARE PHOTOCROMA system that, in addition to increasing the working time for the professional, has an excellent refractive index between the monomers and fillers, and color spreading, perfectly achieving the copy of the dental substrate.

UNICHORMA ACTS... making the dentist's life easier.

DENTISTRY AND PEDIATRIC DENTISTRY

Agility facilitated by rapid restoration due to the bulk increment of up to 5 mm and universal color, thus the service is lighter, calmer and much faster, **especially in relation to children**, without having to leave them for an exhausting time, with their mouths open, in the chair.

PHOTOPOLYMERIZATION SYSTEM

The **SMART CARE FOTOCROMA** polymerization system developed by Smart Dent and incorporated into Atos Unichroma resin, provides professionals with longer working time under the reflector light, especially in relation to its main competitors. The system uses blue light and is compatible with all photopolymerizers on the market that have a wavelength between 400 and 500 nm. Incremental polymerization of 2 mm takes 20 s and bulk polymerization of up to 5 mm takes 30 s.

NANO-HYBRID

Its inorganic composition is composed of high-performance silanized fillers, manufactured in Germany, such as Barium-Aluminum glass in the 0.7 μm range, Silica nanoparticles in the 35 nm range and mainly **Zirconium Nano Silicate** in the 100 nm range, which maximize wear resistance. The uniform distribution of the particles from 35 nm to 0.7 μm contributes to high mechanical properties, excellent handling, high aesthetics and especially shine when polished.

FLUORESCENCE, OPALESCENCE AND RADIOPACITY

Fluorescence is reproduced by special pigments that absorb energy from UV light and emit it through visible light. The Atos Unichroma system is balanced to reproduce this effect, masterfully imitating the fluorescence of natural teeth. The **radiopacity** of the Atos Unichroma system, in addition to being in accordance with ISO 4049, safely provides an accurate diagnosis to the dentist. **Opalescence** is an exclusive characteristic of the enamel in natural teeth that absorbs blue colors and allows orange colors to pass through. For more pronounced opalescence, we recommend using Atos Unichroma with Atos in the color OPL.

STRUCTURING MONOMERS


The organic matrix of this resin is formed by the main structuring monomers in the development of composites in the world, manufactured in the USA and with high molecular weight. They are Bis-Phenol A di-Glycidyl Methacrylate (BisGMA), 1, 12-dodecanediol dimethacrylate (DDDMA), Bis-Phenol A di-Glycidyl Methacrylate ethoxylate (BisEMA), UDMA (Urethane dimethacrylate) and with a slight diluent part of Triethylene glycol dimethacrylate (TEGDMA). These monomers perform important functions, such as reduced polymerization shrinkage factor, color stability, increased service life, thixotropy and ease of handling.

MATERIAL RESISTANCES

The combination of nanometric particles with micrometric particles in the formulation reduces the gap space of the filler particles in the monomeric solution, increasing the percentage of load that, inserted into the structuring monomers, they raise considerably the physical properties and mechanical and resin, providing high resistance, restorations lasting with aesthetic result.

Photomicrograph in Electron Microscopy of Scanning (SEM) of

the resin particle combination. Atos presents uniformity and does not presents discrepant (large) particles. **Magnification of 10,000 times.**
Image kindly provided by FOP-Unicamp

COMPOSITION

Atos Unichroma resin is composed of a matrix with methacrylate monomers, photoinitiator system, silane and stabilizers. Its inorganic filler is composed of Barium-Aluminum glass, Silica nanoparticles and mainly Zirconium Silicate.

UNICHROMATIC

A color for all hues

Light curing time

	Thickness	Time	Color
BULK FILL	5 mm	30 s	UNIVERSAL
INCREMENTAL	2 mm	20 s	

INDICATIONS FOR USE

- Restorations of permanent and deciduous teeth, anterior and posterior, of classes I, II, III*, IV* and V (all classes)**
- Direct composite resin veneers
- Bonding of tooth fragments
- Tooth anchoring (splinting)
- Correction of tooth shape (e.g. conoids)
- Reduction or closure of diastemas
- Non-carious cervical lesions
- Porcelain and/or composite repairs
- Indirect restorations, such as inlays, onlays and veneers

* For restorations with large cavities such as class III and IV (without a back wall to provide color support), use an Opacifying resin (e.g.: Atos resin in the color Opaquer or Dentin A1, A2, etc.) that will serve as a background for the restoration and finish the remaining volume with Atos Unichroma. This procedure is necessary because Atos Unichroma may appear semi-translucent when not surrounded by the walls of the dental cavity. ** See Instructions for Use.

ADVANTAGES

- **Aesthetic dentistry and pediatric dentistry:** Fast restoration with less chair time (no need for color matching).
- **Uncomplicated: no need for layering** opacities such as Enamel, Body and Dentin.
- **Simplicity: no need to know every competitor** between brands, their colors, consistencies, etc.
- **Cost reduction: it is no longer necessary to have stocks of** marks, colors and opacities.
- **Practicality: without the complexity of different color options**
- **and** brands available on the market.

Clinical time saving: Bulk? Incremental? You choose!... Color?

She is unique!

PHYSICOCHEMICAL AND MECHANICAL PROPERTIES

Tested by FACPP and FOP-
UNICAMP universities

POLYMERIZATION DEPTH

The Depth of Polymerization or Depth of Cure measures the efficiency of polymerization and evaluates the clinical longevity of the material. A poor depth of cure causes a risk of pulp aggression due to the toxic characteristics of the unpolymerized monomer, increased wear due to the decrease in the mechanical properties of the resin (mainly its hardness), change in resin color due to incomplete reaction of the material and greater absorption of water and dyes.

The resins were inserted into matrices (5 mm in diameter and 8 mm in depth). The upper portion was protected with a transparent plastic tape in order to press the composite resin and obtain a flat surface. The samples were removed from the matrix after polymerization (Valo, Ultradent, USA) with irradiance maintained at 1200 mW/cm² checked by a radiometer for 40 seconds. The unpolymerized portion of the surface opposite the light source was manually removed using a scalpel blade. The height of the specimen was measured with a digital caliper (Mitutoyo, Japan, with accuracy of 0.01 mm). The value of the polymerization depth was obtained by dividing by two, that is, 50%, the height of the polymerized remnant according to ISO 4049.


Figure 10 Polymerization Depth (mm) and SD in (%) of 5 different composites (n=5). Statistical analysis ANOVA and Tukey test; p<0.05.

Source: Silvestre FA, Kelly K, Moreira ACL, Feitosa VP. Paulo College

Shrike College of Strength (FACPP) and
of Dentistry of
Piracicaba (FOP-Unicamp) 2021.

Atos Unichroma composite resin showed surprising **CURE DEPTH** values, being **SUPERIOR TO ALL** commercial BRANDS of competing resins tested.

DEGREE OF CONVERSION

The degree of conversion of a composite resin is given by how much monomer was converted into polymer, after the polymerization process of the dental composite.

The degree of conversion (DC) of the discs of each adhesive (1 mm x 6 mm) (n = 3) photopolymerized by LED (Valo, Ultradent, South Jordan, USA) with irradiance of 1200 mW/cm² for 40 s, were analyzed by a Fourier transform infrared spectrophotometer (ATR-FTIR, Spectrum Frontier, Perkin-Elmer Corp, Norwalk, USA) in the spectral range of 1650 – 550 cm⁻¹ with a resolution of 4 cm⁻¹ and 16 scans per reading. Moderate pressure (0.034 MPa) was applied to ensure adequate contact between the tooth and the ATR device.


Figure 2: Degree of Conversion (mean and SD in %) of 5 different composites (n=5). Statistical analysis ANOVA and Tukey test; p<0.05.

Source: Silvestre FA, Kelly K, Moreira ACL, Feitosa VP. Paulo College Shrike College of Strength (FACPP) and of Dentistry of Piracicaba (FOP-Unicamp) 2021.

Atos Unichroma is statistically **EQUAL TO ITS MAIN COMPETITORS** of international scope Filtek One (3M) and Tetric N-Ceram (Ivoclar).

FLEXIBLE STRENGTH

Flexural Strength is one of the most important mechanical tests for dental composites as it assesses the fracture resistance of the restorative material as well as

your properties intrinsic (monomeric matrix, loads inorganic and photoinitiators) after photopolymerization, in relation to the different masticatory forces.

The test was carried out in accordance with ISO 4049:2019, using the EMIC-Instron universal testing machine (Model 23-2S, São José dos Pinhais-PR Brazil), with the help of an Odeme split matrix, 5 test specimens were prepared for each resin group with specific dimensions (25x2x2 mm) in the shape of a rod. After 3 photoactivations of 20 s (one on each side and one in the center), with irradiance of 1200 mW/cm², the resins were subjected to a three-point bending test.


Figure 3: Degree of Conversion (mean and SD in %) of 5 different composites (n=5). Statistical analysis ANOVA and Tukey test; $p < 0.05$. **Source:** Silvestre FA, Kelly K, Moreira ACL, Feitosa VP. Faculdade Paulo

Shrike of Strength (FACPP) and College of Dentistry of Piracicaba (FOP-Unicamp) 2021.

Unichroma Acts presents **EXCELLENT FLEXIBLE RESISTANCE**, comparable to well-known products on the market.

COMPRESSION STRENGTH

The mechanical property of compressive strength demonstrates the ability of the composite resin to withstand vertical stresses, since the forces transmitted to the restorations during the chewing act can trigger fractures of the material. The higher the value of compressive strength, the greater the longevity of the composite. The mechanical compression test was performed on the EMIC-Instron universal testing machine (Model 23-2S, São José dos Pinhais-PR Brazil), with a load cell of 50,000 Kgf and speed of 1mm/min. EMIC indicated the total breaking strength of the test specimens. For this purpose, a bluhill version 2.0 program was used (EMIC, São José dos Pinhais - PR, Brazil), which archives the results in a database in the form of graphs and numerical values in MPa.


Figure 4: Degree of Conversion (mean and SD in %) of 5 different composites (n=5). Statistical analysis ANOVA and Tukey test; $p < 0.05$. **Source:** Silvestre FA, Kelly K, Moreira ACL, Feitosa VP. Faculdade Paulo

Shrike of Strength (FACPP) and College of Dentistry of Piracicaba (FOP-Unicamp) 2021.

Atos Unichroma showed balance in its composition, and it was observed that only in this study all results are similar. Atos Unichroma is always **EQUAL OR BETTER** in the other tests.


EXPERT REPORTS

“The tests of **Atos Unichroma** resin showed extremely encouraging results! It worked very well for the **chameleon effect** proposal. A national technology with **excellent performance**. ”

PROFESSOR RENATO VOSS ROSA – CURITIBA (PR)

As with Atos resin, what impressed me most about **Atos Unichroma** was its **excellent consistency**. In addition, its ability to mimic the color of the remaining resin ensures **simplicity for use in Dentistry and Pediatrics**

PROF. DR. SIMONE ALBERTON – PASSO FUNDO (RS)

“**Atos Unichroma** resin has an **impressive ability to mimic the colors** of different substrates. Its handling and clinical application are decisive factors in choosing this resin, without the need for multiple colors. ”

PROF. DR. VICTOR FEITOSA - FACPP PROFESSOR – FORTALEZA (CE)

CLINICAL CASES

RESTORATION WITH ATOS UNICHROMA COMPOSITE RESIN OF FRACTURED CENTRAL INCISOR

Author: Prof. Dr. Fabio Sene - Prof. UEL – Londrina-PR

1. Initial clinical case

2. Application of acid attack

3. Beginning of restoration with Atos
Unichroma

4. Incremental restore

5. Final clinical case - polished and
hydrated

“...simplicity... effectiveness...
great color integration!!!”

PROF. DR. FABIO SENE

ATOS UNICHROMA BULK FILL - SIMPLICITY, EFFICIENCY AND AESTHETICS IN POSTERIOR TEETH

Author: Prof. Dr. Carlos Francci - Prof. from the Department of Biomaterials and Oral Biology at FOUSP

1. Initial case

2. Preparation and selective acid etching

3. Beginning of restoration after application of adhesive

4. Photopolymerization

5. Restored tooth

6. Clinical case completed - final polished and hydrated

PROF. DR. CARLOS FRANCCI

“ Smart Dent surprises by launching a composite that brings together all the latest innovations, Atos Unichroma Bulk Fill, in addition to not needing stratification, it is also unichromatic, that is, it brings together all the simplifications in a single material...
clinical result is unique... congratulations to R&D!!! ”

DIRECT RESTORATION WITH ATOS UNICHROMA IN CAVITY CL II - THE UNION OF SIMPLICITY AND EFFICIENCY.

Author: Prof. Dr. Gustavo Passarelli Petris - Professor in postgraduate and immersion courses at F1 Courses

1. Initial case

2. Initial preparations

3. Installation of the My QuickMat Forte - Polydentia matrix system for shaping the proximal surface.

4. Construction of the proximal box with Atos enamel color A2 simplifying the cavity.

5. Single increment insertion of Atos Unicroma resin.

6. Final result

7. Final result

A Unichroma tem
great sculptability and
versatile features for
various clinical
applications, in
addition to a great and
easy polishing.

PROF. DR. GUSTAVO PASSARELLI PETRIS

SMART DENT

ALSO KNOW

Photo: Prof. Dr. Renato Voss
Rosa

ATOS is a light-curing nano-hybrid composite resin, indicated for restorations of all classes in anterior and posterior teeth using the incremental system.


Photo: Prof. Dr. Maciel
Junior

ATOS ACADEMIC is the colored version of Atos resin, made specifically for teaching techniques dental restoratives.


Photo: Prof. Dr. Idulton
Grabowski

SMART ORTHO is an orthodontic adhesive, developed for direct bonding of metal, ceramic and polycarbonate brackets to the dental surface.

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