UNIK|X

V E N E E R

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Smart Dent - USA

UNIVERSITY OF NORTH CAROLINA
9201 UNIVERSITY CITY BLVD, CHARLOTTE, NC
28223, USA

INTERNATIONAL
TRADE
ADMINISTRATION

UNIK|X

V E N E E R

HIGH PERFORMANCE CEMENTATION FOR CERAMIC
LAMINATES AND AESTHETIC RESTORATIONS

UNIK|X^{LC}
FACET CEMENT

UNIK|X^{LV}
PASTA LOW VISCOSITY

UNIK|X
TRY-IN

Smart Dent - USA
UNIVERSITY OF NORTH CAROLINA
9201 UNIVERSITY CITY BLVD, CHARLOTTE,
NC 28223, USA

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REVOLUTION IN AESTHETIC CEMENTATION

UNIKK Veneer LC (Light-Cured) cement, developed by SmartDent, is a **state-of-the-art photo-cured resin cement**, designed to provide **high adhesion, predictable esthetics and ease of use in the cementation of veneers**, and **UNIKK (LV)** for onlays, inlays and overlays with a maximum thickness of 2 mm.

It stands out from its competitors due to the combination of **innovative technologies that guarantee natural optical behavior, adjusted fluorescence and excellent rheology for passive seating of the restorative piece.**

MAIN DIFFERENTIATORS OF UNIKK VENEER CEMENT

1 - SMART COLOR RANGE - SIMPLIFIED VITA SCALE

While many cements work with multiple complex variations, UNIKK Veneer uses a minimalist and efficient approach **with only 5 colors, uses conventional nomenclatures (BL2, B1, A1, A2 and A3.5) and is translucent**, making it easier to choose the color without the need to "decipher" complex codes.

Easy to choose: If the patient's tooth is A1, simply use A1 cement. No need for confusing codes.

INDICATIONS

- TRANSLUCENT** High opalescence and fluorescence.
- BL2** For whitened teeth, without artificial blue effect.
- B1** The most used in facets and laminates.
- A1 AND A2** For cases of natural color.
- A3.5** For darker substrates.

REVOLUTION IN AESTHETIC CEMENTATION

2 - TECNOLOGIA O.L.D. (OPTIMIZED LIGHT DISPERSION) OPTIMIZED LIGHT DISPERSION


Light scattering within cement is influenced by the size, shape and composition of the material's particles. In UNIKK Veneer, these particles have been tuned to act as a **microscopic prism**, allowing light to **scatter more naturally within the cement** and **pass through to the substrate**.

White light (composed of the red, green and blue spectrums – RGB) enters the cement.

The special particles in the cement act like prisms, opening up the spectrum of light to create a smoother transition between the veneer and the tooth.

A second particle scatters this light in a controlled manner, reducing the risk of shadows or color variations.

The result is a more homogeneous and natural cementation, without that grayish effect that can occur with traditional cements.

3 - FLUORESCENCE IN UNIKK VENEER CEMENT NATURALNESS IN ANY LIGHTING

Fluorescence is the ability of a material to absorb ultraviolet (UV) light and re-emit that light in the visible spectrum, similar to **the behavior of natural tooth enamel**.


The problem: **Many modern ceramics do not fluoresce enough because current ceramic materials were designed to fluoresce under old-fashioned lights (385 nm - Wood's lamps), but modern environments use UV LED lights (405 nm), where the ceramic does not fluoresce as it should.**

Visible fluorescence under modern UV LED (405 nm LED) and UV light (385 nm)

Correction of lack of fluorescence in ultra-thin ceramic veneers

Easier removal of excess with fluorescent microscopy

Predictable and uniform aesthetics in any lighting condition

REVOLUTION IN AESTHETIC CEMENTATION

OPTICAL EFFECT IN CERAMIC LAMINATES: UNDERSTANDING THE INFLUENCE OF THICKNESS AND CEMENT

The choice of cement and the optical interaction between the piece, substrate and cementation material are essential factors for a predictable aesthetic result in restorative dentistry. We will address the main optical effects that occur in ceramic laminates of different thicknesses and how the choice of cement influences the final color of the restoration.

CONTACT LENS EFFECT

(Thickness < 0.5 mm).

Described by Friedman in 1991, the contact lens effect occurs when the piece has minimal thickness and low optical density, allowing the color of the substrate to be visible through it. In these cases, the cement must be translucent, with high fluorescence and opalescence, ensuring the continuity of the passage of light without interfering with the final color.

CHROMATIC EFFECT

(Thickness between 0.5 and 0.7 mm)

When the piece has an intermediate thickness, its color begins to alter the perception of the substrate, creating a final color resulting from the sum of substrate, piece and cement. Here, the choice of cement should be made based on the color of the piece, as its influence is already more significant.

GLASS EFFECT

(Thickness between 0.7 and 0.9 mm)

In this thickness range, there is a high risk of the piece turning gray, a phenomenon known as the glass effect. This problem occurs when the material of the piece, combined with the lack of control in the choice of substrate and cement, results in an undesirable final color. The cement can still interfere with the color, but the main factor becomes the ceramic itself.

LIMITED ROLE OF CEMENT

(Parts Above 0.9mm)

In thicker restorations (over 0.9 mm), the color of the cement has minimal or no influence on the final aesthetics. If the piece already presents an inadequate shade during the test with translucent gel, the ideal solution is to redo it in the laboratory, as no cement will be able to correct this chromatic variation.

VERSATILITY OF UNIKK VENEER CEMENT IN PRACTICE

Smile aesthetics is one of the main challenges in modern dentistry. **This clinical case presents the aesthetic rehabilitation of a patient through the cementation of ceramic laminates, using UNIKK Veneer cement**, a high-tech material with aesthetic predictability.

CASE CHALLENGE

The goal was to correct the shape and contour of the patient's teeth, ensuring a natural and harmonious aesthetic. The main challenge was the differences in thickness of the restorative pieces and choosing the correct cement for each substrate.

Tooth 11 - Class 4 color change and significant darkening;


Teeth 21, 22 and 12 - Color improvement, but need for optical adjustment;

Canine 13 - Small fracture on the incisal edge, requiring an ultra-thin ceramic fragment;

Difference in piece thicknesses - From ultra-thin fragments to conventional preparations requiring a different cementation strategy.

PLANNING AND PREPARATION

The clinical approach followed meticulous planning to define the types of preparation and the choice of cement appropriate for each piece.

A - Ultra-thin fragment (canine 13) - "contact lens" effect for natural transparency;

B - Thin laminates (0.6 to 0.7 mm) - Chromatic effect;

C- Conventional laminate (0.9 mm or more) - Piece with its own color, zero optical effect.

CEMENT SELECTION AND OPTICAL STRATEGY

Cementation was done using UNIKK Veneer, applying different techniques according to the optical needs of each tooth:

- **Canine 13** (Ultra-thin Fragment - Contact Lens Effect)

Translucent Cement - Maintains the passage of light and preserves the original color of the tooth;

- **Thin Laminates** (0.6 to 0.7 mm - Chromatic Effect)

B1 Cement - Harmonizes the color between the piece and the tooth;

- **Tooth 11** (0.9mm Laminate - Zero Effect)

LV (Low Viscosity) Cement - Allows for better light dispersion, ensuring naturalness;

VERSATILITY OF UNIKK VENEER CEMENT IN PRACTICE

CEMENTATION EXECUTION

The application of the cement followed the adhesion and photopolymerization protocols appropriate for each substrate, ensuring safe adhesion and predictable aesthetics.

FINAL RESULT AND CONCLUSION

The final result demonstrates impeccable aesthetic harmony, with a uniform color transition between substrate, cement and ceramic. The use of UNIKK Veneer was essential to ensure:

- **Precise optical correction;**
- **Ideal fluidity in cementation** (P.A.R technology, - Pressure Adaptive Rheology);
- **Homogeneous light dispersion** (O.L.D. technology - Optimized Light Dispersion);
- **Easy removal of excess;**
- **Predictability of the final cementation color.**

Conclusion: UNIKK Veneer cement provided an impeccable finish, correcting chromatic differences and ensuring a natural and harmonious smile. Choosing the correct cement for each type of piece was essential to achieve this level of aesthetic excellence.

CONTACT LENS EFFECT (Thickness < 0.5 mm).

INDICATION OF CEMENTS IN RELATION TO THE THICKNESS OF CERAMIC LAMINATES

UNIKK
TRANSPARENT

CHROMATIC EFFECT (Thickness between 0.5 and 0.7 mm)

UNIKK A1

UNIKK B1

UNIKK BL2

UNIKK A2

UNIKK A3.5

GLASS EFFECT (THICKNESS BETWEEN 0.5 AND 0.7 MM)

UNIQUE A1

UNIQUE B1

UNIQUE BL2

LIMITED ROLE OF CEMENT -UNIKK VENEER LV - POSTERIOR A2 (Parts above 0.9mm < than 2mm)

UNIKK LV B1
(Low Viscosity)

UNIKK VENEER LV - POSTERIOR A2

THE REVOLUTION OF LOW VISCOSITY CEMENTATION

UNIKK Veneer LV Posterior A2 (Low Viscosity) is a photoactivated resin cement, specially developed to ensure passive setting, easy removal of excess and high optical integration. Its advanced technology provides adaptive fluidity under pressure, ensuring precise and predictable cementation, ideal for onlays and overlays with less than 2 mm thickness.

It is the ideal choice for high-performance adhesive cementation. Its intelligent viscosity technology improves the adaptation of the part, prevents cementation failures and ensures a natural and long-lasting result.


1 - BLENDING EFFECT - CHAMELEON EFFECT INTELLIGENT CHROMATIC ADAPTATION

UNIKK Veneer LV has a chameleon effect (Blending Effect), the composition and structure of the cement particles, which were formulated to allow controlled scattering of light, allowing it to **harmonize with the color of the dental structure and the restorative material.** This reduces the need for multiple shades of cement, making the choice simpler and more efficient.

Light entry: The light that falls on the restorative material passes through the cement and interacts with its optical scattering particles.

Chromatic harmonization: The cement particles absorb and reflect light uniformly, avoiding abrupt differences in tone.

Visual integration: The restoration and the tooth become a single optical structure, resulting in a natural and homogeneous appearance.

2 - P.A.R. TECHNOLOGY (PRESSURE ADAPTIVE RHEOLOGY) INTELLIGENT FLUIDITY UNDER PRESSURE

UNIKK Veneer LV is a cement that **combines adaptive fluidity, ease of use and excellent aesthetic integration**, making it the **ideal choice for high-performance adhesive cementation**. Its **intelligent viscosity** technology improves the adaptation of the part, prevents cementation failures and ensures a natural and long-lasting result.

"CHANGE IN CONSISTENCY IS DUE TO THE PRESSURE EXERTED ON THE MATERIAL, HEATING IS NOT REQUIRED"

1 - BEFORE APPLICATION

It has a firm consistency, without dripping, making it easier to position the piece.

APPLICATION OF RESIN CEMENT

2 - DURING APPLICATION

When pressed, the cement becomes fluid and adaptable.

THE NEW TECHNOLOGY OF THIS CEMENT PREVENTS IT FROM RUNNING DURING HANDLING OF THE PIECE.

3 - DURING CEMENTATION

When pressed, the cement becomes fluid and adaptable, allowing for tension-free laying with minimal film thickness.

PART POSITIONING AND PROGRESSIVE MANUAL PRESSURE.

NOTE THE TECHNOLOGY OF UNIKK VENEER LV THAT MODIFIES ITS RHEOLOGY UNDER PRESSURE INCREASING ITS FLOW

4 - AFTER SETTING

It returns to its initial consistency, making it easier to remove excess without effort.

AT THE END OF THE LAYING PROCESS, THE CEMENT RETURNS TO A FIRM CONSISTENCY, ALLOWING THE EXCESS TO BE EASILY REMOVED.

UNIKK VENEER TRY-IN


ABSOLUTE PRECISION IN COLOR SELECTION

Choosing the correct cement color is one of the biggest challenges when cementing veneers and contact lenses. Small color variations can compromise the final aesthetics, making the result unsatisfactory. **To solve this issue, UNIKK Veneer offers one of the most accurate TRY-IN systems on the market, ensuring that the color chosen in the test is identical to the final color after polymerization.**

WHAT IS TRY-IN AND WHY IS IT ESSENTIAL?

TRY-IN is a paste that simulates the color of the cement before its final application. It allows the dentist to test different color options before cementation, **without compromising the piece or the tooth structure.** This step is crucial to ensure **a predictable and natural result, especially in ultra-thin restorations where the color of the cement directly influences the final aesthetics.**

Another innovation of UNIKK Veneer TRY-IN is its compatibility with the Vita color scale, making it easier to choose the right color for each case. While some cements use complex nomenclatures and codes that are not very intuitive.

UNIKK VENEER TRY-IN ABSOLUTE PRECISION IN COLOR SELECTION

COLOR ACCURACY WITH DELTA E < 1 - WHAT DOES THIS MEAN?

When we talk about color accuracy in dentistry, especially in the cementation of veneers and indirect restorations, one of the main challenges is ensuring that the final color of the cement is exactly what the dentist expects.

Delta E (ΔE) is a metric used to measure the difference between two colors. It quantifies how noticeable that difference is to the human eye.

Comparison of Delta E and human perception:

$\Delta E = 0 - 1 \rightarrow$ Imperceptible difference.

$\Delta E = 1 - 2 \rightarrow$ Slightly noticeable difference under rigorous testing conditions.

$\Delta E = 2 - 3 \rightarrow$ Visible difference, but still acceptable.

$\Delta E > 3 \rightarrow$ Easy to notice difference and can compromise aesthetics.

TRY-IN VERSUS FINAL CEMENTATION - PROVEN FIDELITY

Unlike traditional cements, where the color of TRY-IN can change dramatically after polymerization, UNIKK Veneer TRY-IN maintains absolute fidelity between the tested color and the final color. This allows for a safer cementation process, reducing the need for adjustments or rework.

If you have ever had problems with color changes after cementation, UNIKK Veneer solves this problem with an absolutely precise TRY-IN ($\Delta E < 1$). This provides greater safety, predictability and impeccable aesthetics for your patients.

MAIN ADVANTAGES OF USING UNIKK VENEER TRY-IN

MAXIMUM COLOR PREDICTABILITY ($\Delta E < 1$).

EXACT MATCH TO VITA SCALE FOR EASY SELECTION.

REDUCTION OF CLINICAL TIME BY AVOIDING ERRORS IN CEMENT SELECTION.

EASE OF AESTHETIC TESTING

IDEAL FOR ULTRA-THIN VENEERS, ONLAYS AND INLAYS, WHERE THE COLOR OF THE CEMENT INFLUENCES

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